

INNOWALK

In & out

Quick Guide



Pre-check before use

- ✓ Tighten screws and parts if they are loose.
- ✓ Verify the condition of the guide-string.
- ✓ Test the seat and upright function.
- ✓ Contact Made for Movement or its representative if there are any unusual sounds or broken parts..

Transfer the user into the Innowalk

1



Open all straps. Place Innowalk in sitting position. Transfer the user into the Innowalk. It is recommended to use a personal lifter.

2



Place the user well back on the seat.

3



Secure the chest strap.

4



Place the user's leg in the calf bow and attach the guidestring.

NB



For Innowalk Large: Attach the strap for the anti-overstretch on the calf bow.

5



Position the heel towards the back of the foot plate.



Secure with straps.

6 NOTE! Adjust tension of the guide-string.

Innowalk Small



Release tension of the guide-string by using the free-coupling.

Innowalk Medium



Push the red button and release the lock. Adjust the tension by using the handle. Tension increases by turning the handle down to present mark.


Innowalk Large



Turn the handle into position according to the stop-screw. Tension increases by turning the handle to the right.

Note! When adjusted, it is important to keep the same position every time the product is in use.

7



Hold on to the guide-string to control the tension when moving from sitting to standing.

8



Secure the hip belt.

✓ Start the movement and adjust the speed to the user's needs.

✓ Make sure that the user is standing and moving in an optimal position.

Attach necessary accessories

Table



Attach and adjust the table.

Anti-overstretch



Attach the string for the anti-overstretch on the calf bow.

Shoulder straps



Attach shoulder straps.

Handles



Handle for arm motion.

Our highly skilled and experienced consultants advise in optimal transfer techniques.

NOTE!

All people handling the Innowalk should be trained by a Made for Movement representative.

For more information, please see the user manual or contact a Made for Movement representative.


Transfer the user out of the Innowalk

1



Remove the table, shoulder straps, strap for anti-overstretch and **release the hip support**.

2



Hold on to the guide-string to control the tension when moving from standing to sitting.

3



Release the straps over the feet.

4



Release the guide-string and the calfbow.

5



Release the chest support and transfer the user out of the Innowalk. It is recommended to use a personal lifter.



**For more information please
contact us**

 + 47 35 50 51 20

 info@madeformovement.com

 www.madeformovement.com

 **made for
movement**